



Product Code	32246
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

GTIN	00707387790332
------	----------------

## Description

### 3M™ LITTMANN® CLASSIC III STETHOSCOPE

- 3M Code: 5864
- Colour: champagne - burgundy
- Length: 69 cm
- Head Ø: 43/33 mm

The 3M™ Littmann® Classic III™ Stethoscope offers high acoustic sensitivity for exceptional performance when doing general physical assessments. It features dual tunable diaphragms and updated design.

- versatile adult and pediatric two-sides chestpiece
- dual single-piece tunable diaphragms (pressure sensitive)
- small-side diaphragm prevents debris in open bell
- high acoustic sensitivity
- snap tight soft-sealing eartips
- next generation longer-life tubing
- small side converts to open bell