



Product Code	32250
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

EAN/UPC	4250296163903
UMDNS	13-755

Description

ERKA FINESSE STETHOSCOPES

FINESSE LIGHT ADULT - head Ø 44.5/33.5 mm
ERKA code: 520 000 10
Finesse light - black
tubing length - 55 cm

With the reversible double chest-piece you'll be prepared for all eventualities in everyday medical practice. The membrane side guarantees optimal auscultation for the entire frequency spectrum and the open bell side is specifically designed for low frequencies.

FEATURES

Pre-tensioned ERKA convex membrane provides excellent resonance and sound.
Superior polish, chromium-plated, nonporous surface and brass chest-piece ensure perfect surface hygiene.
Ergonomic design and super soft tips provide excellent acoustic seal against external noise interference and maximize user's comfort.
Standard 15° inclined binaurals can be individually adapted to the specific ear canal of the user.

All Erka stethoscopes are supplied with carrying bag

Made in Germany.