



Product Code	32258
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	4250296148672
UMDNS	13-755

## Description

### ERKA FINESSE STETHOSCOPES

FINESSE PEDIATRIC - head Ø 31/26.5 mm  
ERKA code: 549 000 35  
Finesse pediatric - pink  
tubing length - 55 cm

With the reversible double chest-piece you'll be prepared for all eventualities in everyday medical practice. The membrane side guarantees optimal auscultation for the entire frequency spectrum and the open bell side is specifically designed for low frequencies.

#### FEATURES

Pre-tensioned ERKA convex membrane provides excellent resonance and sound.  
Non-chill ring ensures highest patient comfort.  
Superior polish, chromium-plated, nonporous surface and brass chest-piece ensure perfect surface hygiene.  
Two internal channels run parallel in the stethoscope tubing preventing any noise caused by friction of two separate tubes.  
Through the separately integrated spring in the tube, there is no acoustic loss from the chest-piece to the ear.  
Ergonomic design and super soft tips provide excellent acoustic seal against external noise interference and maximize user's comfort.  
Standard 15° inclined binaurals can be individually adapted to the specific ear canal of the user.

All Erka stethoscopes are supplied with carrying bag

Made in Germany.