



Product Code	32270
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	4250296148245
UMDNS	13-755

## Description

### ERKA FINESSE2 STETHOSCOPES

FINESSE2 PEDIATRIC - head Ø 30/21.3 mm  
ERKA code: 536 000 00  
Finesse2 pediatric - black  
tubing length - 55 cm

The reversible, combined chest-piece has two membrane sides, each making the full frequency range perfectly audible. The large membrane is suitable for auscultation of adults and the small one for children or slim patients.

### FEATURES

Pre-tensioned ERKA convex membrane provides excellent resonance and sound.  
Non-chill ring ensures highest patient comfort.  
Superior polish, chromium-plated, nonporous surface and brass chest-piece ensure perfect surface hygiene.  
Two internal channels run parallel in the stethoscope tubing preventing any noise caused by friction of two separate tubes.  
Through the separately integrated spring in the tube, there is no acoustic loss from the chest-piece to the ear.  
Ergonomic design and super soft tips provide excellent acoustic seal against external noise interference and maximize user's comfort.  
Standard 15° inclined binaurals can be individually adapted to the specific ear canal of the user.

All Erka stethoscopes are supplied with carrying bag

Made in Germany.