



| | |
|---------------|----------------|
| Product Code | 32290 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------|-------|
| CND | C9005 |
| GMDN | 13755 |

Description

3M™ LITTMANN® CLASSIC III STETHOSCOPE

- 3M Code: 5620
- Colour: black
- Length: 69 cm
- Head Ø: 43/33 mm

The 3M™ Littmann® Classic III™ Stethoscope offers high acoustic sensitivity for exceptional performance when doing general physical assessments. It features dual tunable diaphragms and updated design.

- versatile adult and pediatric two-sides chestpiece
- dual single-piece tunable diaphragms (pressure sensitive)
- small-side diaphragm prevents debris in open bell
- high acoustic sensitivity
- snap tight soft-sealing eartips
- next generation longer-life tubing
- small side converts to open bell