

LITTMANN CLASSIC III - 5910C - champagne rose - stainless finish



Product Code	32321
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

GTIN 00707387814151

Description

3M™ LITTMANN® CLASSIC III STETHOSCOPE SPECIAL EDITION

- 3M Code: 5910C
- Colour: champagne rose - stainless finish
- Length: 69 cm
- Head Ø: 43/33 mm

The 3M™ Littmann® Classic III™ Stethoscope offers high acoustic sensitivity for exceptional performance when doing general physical assessments. It features dual tunable diaphragms and updated design.

- versatile adult and pediatric two-sides chestpiece
- dual single-piece tunable diaphragms (pressure sensitive)
- small-side diaphragm prevents debris in open bell
- high acoustic sensitivity
- snap tight soft-sealing eartips
- next generation longer-life tubing
- small side converts to open bell