

RIESTER DUPLEX 2.0 ALUMINIUM STETHOSCOPE - adult - all black



Product Code	32329
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2741198
CND	C9005
EAN/UPC	4045396195424

Description

RIESTER DUPLEX® 2.0 ALUMINIUM STETHOSCOPES - LATEX FREE

- Riemer code: 4201-01
- Adult
- Colour: all black
- Membrane: Ø 32/44 mm
- Length: 74 cm
- Weight: 97 g

Duplex® 2.0 stethoscopes provide outstanding acoustic performance. Special design tube for better auscultation of adults, children/babies or new-borns.

- precision double chest-piece made of stainless steel with innovative membrane support
- anatomically shaped binaural with internal multiple spring for adjusting the contact pressure
- extra-soft, replaceable and pivotable ear-tips for perfect sealing at the ear canal and enhanced wear comfort
- innovative acoustic system with extraordinary auscultation for all frequency ranges
- special membrane for precise acoustics with novel non-chill rims for improved adaptation on the skin

Supplied with name plate, replacement ear-tips and membrane.

Made in Germany.