

RIESTER DUPLEX 2.0 S/S STETHOSCOPE - adult - red



Product Code	32336
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1158487
CND	C9005
EAN/UPC	4045396195462
GMDN	13755

Description

RIESTER DUPLEX® 2.0 STAINLESS STEEL STETHOSCOPE - ADULT

- Riemer code: 4210-04
- Colour: red
- Membrane: Ø 32/44 mm
- Length: 74 cm
- Weight: 151 g
- Latex free

Duplex® 2.0 stethoscopes are distinguished by outstanding acoustic performance. They have been completely new constructed to meet the challenges of daily diagnostic practice. Special design tube for better auscultation of adults, children/babies or new-borns.

- precision double chest-piece made of stainless steel with innovative membrane support
- anatomically shaped binaural with internal multiple spring for adjusting the contact pressure
- extra-soft, replaceable and pivotable ear-tips for perfect sealing at the ear canal and enhanced wear comfort
- innovative acoustic system with extraordinary auscultation for all frequency ranges
- special membrane for precise acoustics with novel non-chill rims for improved adaptation on the skin

Supplied with name plate, replacement ear-tips and membrane.

Made in Germany.