

# RIESTER DUPLEX 2.0 S/S STETHOSCOPE - pediatric - white



Product Code	32338
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1158546
CND	C9005
EAN/UPC	8023279323382
GMDN	13755

## Description

### RIESTER DUPLEX® 2.0 STAINLESS STEEL STETHOSCOPE - PEDIATRIC

- Riester code: 4220-02
- Colour: white
- Membrane: Ø 32/23.5 mm
- Length: 69 cm
- Weight: 124 g
- Latex free

Duplex® 2.0 stethoscopes are distinguished by outstanding acoustic performance. They have been completely new constructed to meet the challenges of daily diagnostic practice. Special design tube for better auscultation of adults, children/babies or new-borns.

- precision double chest-piece made of stainless steel with innovative membrane support
- anatomically shaped binaural with internal multiple spring for adjusting the contact pressure
- extra-soft, replaceable and pivotable ear-tips for perfect sealing at the ear canal and enhanced wear comfort
- innovative acoustic system with extraordinary auscultation for all frequency ranges
- special membrane for precise acoustics with novel non-chill rims for improved adaptation on the skin

Supplied with name plate, replacement ear-tips and membrane.

Made in Germany.