

ERKAPHON FLAT CHEST-PIECE - aluminium - black 544 000 10



Product Code	32344
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1146891
CND	C9005
EAN/UPC	4250296129640
GMDN	13755

Description

ERKAPHON STETHOSCOPES

ERKAPHON FLAT CHESTPIECE - aluminium - black
ERKA code: 544 000 10

High quality precision stethoscopes for excellent acoustics.
All stethoscopes are equipped with 56 cm long Y tube, aluminium binaural and chest-piece, membrane Ø 44 mm.

Made in Germany.