



Product Code	32350
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	1146910
CND	C9005
EAN/UPC	4250296129466
GMDN	13755

Description

ERKAPHON STETHOSCOPIES

ERKAPHON FLAT CHESTPIECE - chromed - black

ERKA code: 543 000 10

High quality precision stethoscopes for excellent acoustics.

All stethoscopes are equipped with 56 cm long Y tube, brass chrome-plated binaural and chestpiece, membrane Ø 44 mm.

Made in Germany.