

LITTMANN "LIGHTWEIGHT II S.E." - 2450 - black



Product Code	32390
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

GTIN	00707387789145
GMDN	13755

Description

3M™ LITTMANN® LIGHTWEIGHT II S.E. STETHOSCOPE

- 3M Code: 2450
 - Colour: black
 - Head: Ø 44 mm
 - Length: 71 cm
-
- Lightweight: only 118 g
 - Headset made of aerospace alloy
 - Tunable diaphragm for easy low and high frequency monitoring