

# LITTMANN "CLASSIC II INFANT" - 2114 - black



|               |                |
|---------------|----------------|
| Product Code  | 32410          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|      |                |
|------|----------------|
| GTIN | 00707387788759 |
| GMDN | 13755          |

## Description

3M™ LITTMANN® CLASSIC II STETHOSCOPE - INFANT

- 3M Code: 2114
- Colour: black
- Lenght: 71 cm
- Head Ø: 27/19 mm