

SHORT OBSTETRIC STETHOSCOPE - wood



Product Code	32501
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279325010
PARAF	903143473

Description

SHORT OBSTETRIC STETHOSCOPE - beech-wood

These stethoscopes are hand-craft manufactured, made of beech-wood accurately selected and seasoned for at least 3 years. They have a perfect finish with a special transparent varnish.
Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, RO, HU, GR, Arabic.

Length 13 cm

Made in Germany.