



Product Code	32506
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512561
CND	C9005
EAN/UPC	8023279325065
PARAF	987473307
GTIN	08023279325065
GMDN	13755

## Description

WAN PLUS double head pediatric - LATEX FREE - black  
 with inner spring binaural, superior tubing quality and 2 pairs of soft colour coordinated ear tips in different sizes.

Features a chrome-plated zinc-alloy chestpiece, a screwon plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring.

Assure optimum acoustic performances.

Tubing length: 63 cm.

Pediatric stethoscope: diaphragm Ø 35 mm, bell Ø 30 mm.

Supplied in a cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.