

WAN NEONATAL STETHOSCOPE - grey



Product Code	32509
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2639347
CND	C9005
EAN/UPC	8023279325096
PARAF	903176042
GTIN	08023279325096
GMDN	13755

Description

WAN double head neonatal - LATEX FREE - grey

Features a chrome-plated zinc-alloy chestpiece, a screw-on plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring.
Assure optimum acoustic performances.

Chestpiece: Ø 22 mm, bell Ø 19 mm.
Tubing length: 63 cm.
Supplied in a cardboard box.
Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.