

# WAN PEDIATRIC DUAL HEAD STETHO - Y red



Product Code	32515
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2529115
CND	C9005
EAN/UPC	8023279325157
PARAF	987473283
GTIN	08023279325157
GMDN	13755

## Description

WAN double head pediatric - LATEX FREE - red

Features a chrome-plated zinc-alloy chestpiece, a screw-on plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring.  
Assure optimum acoustic performances.

Chestpiece: Ø 35 mm, bell Ø 30 mm.

Tubing length: 63 cm.

Supplied in a cardboard box.

Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.