



Product Code	32516
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512572
CND	C9005
EAN/UPC	8023279325164
PARAF	903176865
GTIN	08023279325164
GMDN	13755

## Description

CLASSIC Double Head S/S Stethoscope - pediatric - black

High quality stethoscope with stainless steel chestpiece, stainless steel double inner spring binaural and durable, elegantly moulded PVC tubing for unsurpassed acoustic performance.

The colour coordinated flexible 55 cm (22") PVC "Y" tubing, non-chill ring and diaphragm retaining ring provide excellent patient comfort, while adding attractiveness.

It comes in a foam lined box with 2 extra sets of soft eartips.

Chestpiece: Ø 45 mm (adult), Ø 35 mm (pediatric).

Manual and box in: GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, RO, BG, Arabic.