



Product Code	32524
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512551
CND	C9005
EAN/UPC	8023279325249
GTIN	08023279325249
GMDN	13755

## Description

### LINUX STETHOSCOPE - black

Single head stethoscope with lightweight zinc alloy chestpiece plated with an attractive satin finish.  
 The special design assures better performances by pushing the back of the head.  
 Non-chill diaphragm retaining rim for patient comfort.  
 Stainless steel binaurals fixed at 15° angle for greater comfort and acoustic seal.  
 Soft and flexible PVC tubing.  
 Supplied in deluxe box with one pair of plastic eartips and PVC eartips.  
 Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.

Chestpiece Ø: 44 mm  
 Tubing: length 62 cm  
 int. Ø 4 mm  
 ext. Ø 8 mm