



Product Code	32526
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2639355
CND	C9005
EAN/UPC	8023279325263
PARAF	903177057
GMDN	13755

Description

REGALITE - black

Top quality single head stethoscope.
Chromed plated chestpiece with "SUPERFLEX" inner spring binaural and 1 piece moulded, PVC tubing.

Supplied in deluxe box, with spare soft eartips and one spare diaphragm.
Manual and box in 7 languages: GB, FR, IT, ES, PT, GR, Arabic.