

CLASSIC DUAL HEAD STETHO - Y light blue



Product Code	32533
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512572
CND	C9005
EAN/UPC	8023279325331
PARAF	987471669
GTIN	08023279325331
GMDN	13755

Description

CLASSIC Double Head S/S Stethoscope - ADULT - light blue

High quality stethoscope with stainless steel chestpiece, stainless steel double inner spring binaural and durable, elegantly moulded PVC tubing for unsurpassed acoustic performance.

The colour coordinated flexible 55 cm (22") PVC "Y" tubing, non-chill ring and diaphragm retaining ring provide excellent patient comfort, while adding attractiveness.

It comes in a foam lined box with 2 extra sets of soft eartips.

Chestpiece: Ø 45 mm (adult), Ø 35 mm (pediatric).

Manual and box in: GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, RO, BG, Arabic.