

CLASSIC DUAL HEAD STETHO - Y burgundy



Product Code	32535
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512572
CND	C9005
EAN/UPC	8023279325355
PARAF	903177273
GTIN	08023279325355
GMDN	13755

Description

CLASSIC Double Head S/S Stethoscope - ADULT - burgundy

High quality stethoscope with stainless steel chestpiece, stainless steel double inner spring binaural and durable, elegantly moulded PVC tubing for unsurpassed acoustic performance.

The colour coordinated flexible 55 cm (22") PVC "Y" tubing, non-chill ring and diaphragm retaining ring provide excellent patient comfort, while adding attractiveness.

It comes in a foam lined box with 2 extra sets of soft eartips.

Chestpiece: Ø 45 mm (adult), Ø 35 mm (pediatric).

Manual and box in: GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, RO, BG, Arabic.