

CLASSIC DUAL HEAD STETHO - Y royal blue



Product Code	32540
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2512572
CND	C9005
EAN/UPC	8023279325409
PARAF	903177436
GTIN	08023279325409
GMDN	13755

Description

CLASSIC Double Head S/S Stethoscope - ADULT - royal blue

High quality stethoscope with stainless steel chestpiece, stainless steel double inner spring binaural and durable, elegantly moulded PVC tubing for unsurpassed acoustic performance.

The colour coordinated flexible 55 cm (22") PVC "Y" tubing, non-chill ring and diaphragm retaining ring provide excellent patient comfort, while adding attractiveness.

It comes in a foam lined box with 2 extra sets of soft eartips.

Chestpiece: Ø 45 mm (adult), Ø 35 mm (pediatric).

Manual and box in: GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, RO, BG, Arabic.