



Product Code	32545
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279325454
GTIN	08023279325454

## Description

REGALITE VET - black

Top quality single head stethoscope.  
Chromed plated chestpiece with "SUPERFLEX" inner spring binaural and 1 piece moulded PVC, tubing 90 cm long.

Supplied in deluxe box, with spare soft eartips and one spare diaphragm.

## Technical Specifications

- 90 cm long Y tube