



Product Code	32550
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2529115
CND	C9005
EAN/UPC	8023279325508
PARAF	987471632
GTIN	08023279325508
GMDN	13755

Description

Cardiology classic stethoscope - blue

Cardiology Classic offers acoustical excellence in a traditional design with superior durability.

Its innovative features include:

- stainless steel chestpiece (Ø 33-44 mm) with extra deep cone-shaped bell for a superb low frequency response
- stainless steel double inner spring binaural
- ultra-sensitive diaphragm for better amplification
- durable, elegantly moulded 55.5 cm (22") PVC tubing (int. Ø 6.5 mm, ext. Ø 12 mm)

Manual and box in 11 languages: GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, Arabic.

Colour coordinated non-chill ring and diaphragm retaining ring, PVC tubing.
It comes in a foam lined box with 2 extra sets of PVC eartips.