



Product Code	32559
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2512789
CND	C9005
EAN/UPC	8023279325591
PARAF	926058734
GTIN	08023279325591
GMDN	13755

Description

TRAD STETHOSCOPES - LATEX FREE

- TRAD single head - pediatric - light blue

TRAD STETHOSCOPES: adult, pediatric

Traditional aluminium stethoscopes, lightweight and reliable.

Chestpiece diaphragm Ø 43.5 mm (Ø 35 for pediatric model), bell Ø 30 mm (Ø 26.5 for pediatric model) external sizes.

Brass binaural, non chill ring.

Tubing: length 62 cm, int. Ø 4 mm, ext. Ø 8 mm.

Supplied in cardboard box. Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.