



| | |
|---------------|----------------|
| Product Code | 32568 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2529115 |
| CND | C9005 |
| EAN/UPC | 8023279325683 |
| PARAF | 903178123 |
| GTIN | 08023279325683 |
| GMDN | 13755 |

Description

TRAD STETHOSCOPES - LATEX FREE

- TRAD double head - dark green

TRAD STETHOSCOPES: adult, pediatric

Traditional aluminium stethoscopes, lightweight and reliable.

Chestpiece diaphragm Ø 43.5 mm (Ø 35 for pediatric model), bell Ø 30 mm (Ø 26.5 for pediatric model) external sizes.

Brass binaural, non chill ring.

Tubing: length 62 cm, int. Ø 4 mm, ext. Ø 8 mm.

Supplied in cardboard box. Manual and box in 11 languages: GB, FR, IT, ES, DE, PT, GR, PL, HU, SE, Arabic.