

WAN STETHOSCOPE - Y blue



Product Code	32570
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2639355
CND	C9005
EAN/UPC	8023279325706
PARAF	903178248
GTIN	08023279325706
GMDN	13755

Description

WAN ADULT STETHOSCOPES - latex free

WAN single head adult - royal blue

Features a chrome-plated zinc-alloy chestpiece, a screw-on plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring.

Assure optimum acoustic performances.

Chestpiece: Ø 45 mm

Tubing: length 63 cm, int. Ø 4.2 mm, ext. Ø 8 mm.

Supplied in a cardboard box.

Multilingual box and manual: GB, FR, IT, ES, DE, PT, HR, GR, BG, RO, PL, HU, SE, FI, Arabic.