



| | |
|---------------|----------------|
| Product Code | 32576 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2513110 |
| CND | C9005 |
| EAN/UPC | 8023279325768 |
| PARAF | 903178325 |
| GTIN | 08023279325768 |
| GMDN | 13755 |

Description

WAN ADULT STETHOSCOPES - latex free

WAN double head adult - burgundy

Features a chrome-plated zinc-alloy chestpiece, a screw-on plastic retaining ring, outer spring binaural and colour matched moulded PVC tubing with non-chill ring.

Assure optimum acoustic performances.

Chestpiece: Ø 45 mm, bell Ø 30 mm.

Tubing: length 63 cm, int. Ø 4.2 mm, ext. Ø 8 mm.

Supplied in a cardboard box.

Multilingual box and manual: GB, FR, IT, ES, DE, PT, HR, GR, BG, RO, PL, HU, SE, FI, Arabic.