



Product Code	32584
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2513117
CND	C9005
EAN/UPC	8023279325843
PARAF	903178539
GTIN	08023279325843
GMDN	13755

## Description

Gima Sprague Rappaport type double-head, double-tube, latex free stethoscope at the price of a standard stethoscope.

- 5 in 1, five way convertible type stethoscope

This convertible stetho, zinc alloy treated with chrome plated chestpiece drum fitting with different sized chestpieces bells, provides versatile acoustics for infants and adults for the detection of a wide range of heart and lung sounds. Chestpiece with small and large diaphragm.

- Double tubing for individual sound transmission

Two separate 56 cm (22") thick latex free wall tubing provide individual sound transmission from chestpiece to both sides of the ear.

- Complete in cardboard box with accessory kit (32590)

It is supplied with an accessory pouch containing: 3 interchangeable plastic black non-chill bells, adult and pediatric spare diaphragm, one pair clear neoprene eartips, one pair black mushroom eartips.

Multilingual manual (GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, RO, NL, LV, SK, Arabic).