



|               |                |
|---------------|----------------|
| Product Code  | 32585          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2513117        |
| CND        | C9005          |
| EAN/UPC    | 8023279325850  |
| PARAF      | 987472267      |
| GTIN       | 08023279325850 |
| GMDN       | 13755          |

## Description

Gima Sprague Rappaport type double-head, double-tube, latex free stethoscope at the price of a standard stethoscope.

- 5 in 1, five way convertible type stethoscope

This convertible stetho, zinc alloy treated with chrome plated chestpiece drum fitting with different sized chestpieces bells, provides versatile acoustics for infants and adults for the detection of a wide range of heart and lung sounds. Chestpiece with small and large diaphragm.

- Double tubing for individual sound transmission

Two separate 56 cm (22") thick latex free wall tubing provide individual sound transmission from chestpiece to both sides of the ear.

- Complete in cardboard box with accessory kit (32590)

It is supplied with an accessory pouch containing: 3 interchangeable plastic black non-chill bells, adult and pediatric spare diaphragm, one pair clear neoprene eartips, one pair black mushroom eartips.

Multilingual manual (GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, RO, NL, LV, SK, Arabic).