

JOTARAP DUAL HEAD STETHOSCOPE - Y transparent pink



Product Code	32587
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2513117
CND	C9005
EAN/UPC	8023279325874
PARAF	987472293
GTIN	08023279325874
GMDN	13755

Description

Gima Sprague Rappaport type double-head, double-tube, latex free stethoscope at the price of a standard stethoscope.

- 5 in 1, five way convertible type stethoscope

This convertible stetho, zinc alloy treaded with chrome plated chestpiece drum fitting with different sized chestpieces bells, provides versatile acoustics for infants and adults for the detection of a wide range of heart and lung sounds. Chestpiece with small and large diaphragm.

- Double tubing for individual sound transmission

Two separate 56 cm (22") thick latex free wall tubing provide individual sound transmission from chestpiece to both sides of the ear.

- Complete in cardboard box with accessory kit (32590)

It is supplied with an accessory pouch containing: 3 interchangeable plastic black non-chill bells, adult and pediatric spare diaphragm, one pair clear neoprene eartips, one pair black mushroom eartips.

Multilingual manual (GB, FR, DE, IT, ES, PT, GR, PL, HU, SE, RO, NL, LV, SK, Arabic).