



Product Code	32594
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1783045
CND	Z12040199
GTIN	00850010298009

Description

3M EKO CORE DIGITAL ATTACHMENT - 8481

Eko CORE allows to transform Littmann Cardiology IV from analog into digital stethoscopes.

Single button press digitally amplifies heart, lung and other body sounds up to 40x (at peak frequency vs. analog) adding automated murmur detection (88% sensitivity and specificity), active noise cancellation and wireless listening. Compatible with most major brands and vet stheto.

When ON, Eko CORE turns stethoscope into a smart detection tools, when OFF, stethoscope still works as an analog Cardiology IV.

Supplied with rechargeable lithium battery (8 hour autonomy) and micro USB charging.