

# MEDIKALL CLEAN PROOF S HIGIENIC COVER FOR STETHOSCOPES



Product Code	32595
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	4250296185370

## Description

MEDIKALL CLEAN PROOF S - box of 500

Medikall clean Proof S is proven to be an effective disposable barrier for stethoscopes against germs. It doesn't decrease the quality of auscultation. Easy to use, efficient through significant time saving and first of all: safe for protecting patients and physicians. Multilingual box: GB, FR, IT, DE, ES, PT.

**PROTECTION GUARANTEED**  
By the effective and safe germ barrier

**TIME SAVING**  
Significant time saving compared to conventional disinfection with exposure time

**QUALITY IMPROVEMENT**  
No loss in the quality of auscultation

**SAFETY**  
It avoids risk of cross contamination

Made in Germany.