

# RIGID EARTIPS SCREW TYPE - white for Classic/Wan/Yton



Product Code	32658
--------------	-------

Unit of sale	box of 10 pcs
--------------	------------------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279326581
---------	---------------

PARAF	987469943
-------	-----------

## Description

Rigid eartips screw-type - white  
spare for Classic/Wan/Yton/Duca/Regalite