

CLASSIC CASE for stethoscope - blue



Product Code	32665
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

EAN/UPC	8023279326659
---------	---------------

Description

CLASSIC CASE FOR STETHOSCOPES - blue

High quality protective hard case for classic style, single tube stethoscopes, available in 3 colours. Its design holds the stethoscope secured in place by a preformed interior, keeping the sensitive diaphragm and tubing protected. Lightweight EVA construction and strong protective outer shell covered in 1680D nylon. Pre-formed moulded interior can also be removed for flexible usage. Sealed dustproof and waterproof zipper with rubber easy-grip zipper pulls for easy opening and closing. Interior mesh pocket for additional accessory storage. Carry handle for easy transportation.

Outside size 28x16x5 cm; inside size 26.5x14.5x4 cm.
Weight: 220 g

Preformed eartips holders
Removable EVA tray with multispandex fabric
Mesh pocket with elastic, 14x19 cm.