

CLASSIC CASE for stethoscope - turquoise



Product Code	32669
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279326697

Description

CLASSIC CASE FOR STETHOSCOPES - turquoise

High quality protective hard case for classic style, single tube stethoscopes, available in 3 colours. Its design holds the stethoscope secured in place by a preformed interior, keeping the sensitive diaphragm and tubing protected. Lightweight EVA construction and strong protective outer shell covered in 1680D nylon. Pre-formed moulded interior can also be removed for flexible usage. Sealed dustproof and waterproof zipper with rubber easy-grip zipper pulls for easy opening and closing. Interior mesh pocket for additional accessory storage. Carry handle for easy transportation.

Outside size 28x16x5 cm; inside size 26.5x14.5x4 cm.
Weight: 220 g

Preformed eartips holders
Removable EVA tray with multispandex fabric
Mesh pocket with elastic, 14x19 cm.