



|              |       |
|--------------|-------|
| Product Code | 32670 |
|--------------|-------|

|              |      |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

|               |   |
|---------------|---|
| Minimum order | 1 |
|---------------|---|

|      |                   |
|------|-------------------|
| Type | No medical device |
|------|-------------------|

|       |  |
|-------|--|
| Class |  |
|-------|--|

|         |               |
|---------|---------------|
| EAN/UPC | 8023279326703 |
|---------|---------------|

## Description

### CLASSIC CASE FOR STETHOSCOPES - pink

High quality protective hard case for classic style, single tube stethoscopes, available in 3 colours. Its design holds the stethoscope secured in place by a preformed interior, keeping the sensitive diaphragm and tubing protected. Lightweight EVA construction and strong protective outer shell covered in 1680D nylon. Pre-formed moulded interior can also be removed for flexible usage. Sealed dustproof and waterproof zipper with rubber easy-grip zipper pulls for easy opening and closing. Interior mesh pocket for additional accessory storage. Carry handle for easy transportation.

Outside size 28x16x5 cm; inside size 26.5x14.5x4 cm.

Weight: 220 g

Preformed eartips holders

Removable EVA tray with multispandex fabric

Mesh pocket with elastic, 14x19 cm.