



Product Code	32671
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1147544
CND	C9006
EAN/UPC	4250296151542
UMDNS	16-156
PARAF	987475023

Description

VARIO - desk - adult cuff - 285.20.481

Superior quality product for use in consulting rooms or hospital.

The Vario model consists of a high precision manometer in a high quality ABS housing with an integrated fastening spigot for individual variable swivelling, a big cuff basket in white powder coated sheet steel, a bulb holder and spiral tubing with bulb and air release valve. This high quality chromium-plated valve makes sensitive adjustments easy.

The 150 mm scale is equipped with a mirrored exterior ring and non-reflective manometer glass.

Made in Germany

Technical Specifications

Graduation: 0-300 mmHg ± 3 mmHg

Weight with basket: 1,350 g

Height of trolley: 830-1,360 mm (adjustable)

Tubing system \varnothing : 18-25 mm

Swivel range: up to 60 deg°