



Product Code	32672
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	178522
CND	C9006
EAN/UPC	8023279326727
UMDNS	16-156
PARAF	987475011

Description

VARIO - wall - adult cuff - 288.20.481+295.00.003

Superior quality product for use in consulting rooms or hospital.

The Vario model consists of a high precision manometer in a high quality ABS housing with an integrated fastening spigot for individual variable swivelling, a big cuff basket in white powder coated sheet steel, a bulb holder and spiral tubing with bulb and air release valve. This high quality chromium-plated valve makes sensitive adjustments easy.

The 150 mm scale is equipped with a mirrored exterior ring and non-reflective manometer glass.

Made in Germany

Technical Specifications

Height with basket: 260 mm
 Width with basket: 180 mm
 Depth with basket: 120 mm
 Graduation: 0-300 mmHg ± 3 mmHg
 Weight with basket: 1,350 g
 Height of trolley: 830-1,360 mm (adjustable)
 Tubing system \varnothing : 18-25 mm
 Swivel range: up to 60 deg°