

# LITTMANN CARDIOLOGY IV - 6239 - plum - high polish rainbow finish - violet stem



Product Code	32674
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

PARAF 987474552

## Description

### 3M™ LITTMANN® CARDIOLOGY IV STETHOSCOPE - SPECIAL FINISHES

- Length: 69 cm
- Head Ø: 43/33 mm
- 3M Code: 6239
- Colour: plum - high polish rainbow finish - violet stem

The 3M™ Littmann® Cardiology IV helps you to catch more of the subtle critical sounds, thanks to the improved acoustics and the new design, more versatile than ever.

- better audibility of high-frequency sounds
- tunable diaphragms for hearing both high and low frequency sounds without turning over the chestpiece (pressure sensitive)
- adult and paediatric sides for more versatility
- single-piece diaphragm, easier to attach and clean, because the surface is smooth without crevices
- next generation longer-life tubing
- snap tight soft-sealing eartips.