

# INFUSION CUFF - 500 ml



Product Code	32680
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2732954
CND	B99
EAN/UPC	8023279326802
PARAF	903170470
GTIN	08023279326802
GMDN	13100

## Description

INFUSION CUFF - 500 ml

Infusion cuff with quality manometer up to 300 mmHg.  
Latex bladder and bulb.

Supplied with carrying bag and multilingual instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.