



Product Code	32684
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396138780

Description

RIESTER METPAK® INFUSION CUFF - 500 ml

Pressure cuff for infusion with high quality chrome-plated manometer Ø 49 mm up to 300 mmHg that allows accurate pressure control. Latex free. Sturdy blue cotton cuff washable at 60° C.

Clear mesh front to allow visualization of fluid level

Made in Germany.