



| | |
|---------------|----------------|
| Product Code | 32686 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I M |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

Description

SANAPHON® SPHYGMO with integrated stethoscope

A classic palm-style sphygmomanometer for precise self-measurement including stethoscope. Sanaphon® guarantees exact measurement values even for patients with heart rhythm disturbances. Maximum error tolerance ± 3 mmHg. Supplied in zipper bag.

Sturdy and light aluminium casing, corrosion protected
Spoon-shaped handle made of stainless-steel

Made in Germany.