



Product Code	32686
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

GTIN	04045396002012
------	----------------

Description

SANAPHON® SPHYGMO with integrated stethoscope

A classic palm-style sphygmomanometer for precise self-measurement including stethoscope. Sanaphon® guarantees exact measurement values even for patients with heart rhythm disturbances. Maximum error tolerance ± 3 mmHg. Supplied in zipper bag.

Sturdy and light aluminium casing, corrosion protected
Spoon-shaped handle made of stainless-steel

Made in Germany.