



Product Code	32687
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1924417
CND	C9006
EAN/UPC	4045396178564
PARAF	987474917

Description

E-MEGA® SPHYGMO - white - adult

Disinfectable, washable and latex free one-piece cuff.

- fast inflation of the cuff due to ergonomic spoon design
- lifelong precision due to almost non-aging precision movement, with specially tempered copper-beryllium diaphragm for pressure-loading capacity up to 600 mmHg
- maximum error tolerance ± 3 mmHg
- comfort - air release valve, wearfree, with fine adjustment

Made in Germany.

- High-polished manometer surface, easy to disinfect
- Large dial \varnothing 63 mm
- Supplied in zipper bag
- Shock protecting polycarbonate case