

# KOBE PALM ANEROID SPHYGMOMANOMETER



Product Code	32690
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	2912377
CND	C9006
EAN/UPC	6970384251852
UMDNS	16156
PARAF	903170710

## Description

### KOBE PALM SPHYGMOMANOMETER

Economic, lightweight, palm type sphygmomanometer with reliable large dial manometer for easy reading.  
Scale 0-300 mmHg, graduated every 2 mmHg.  
Adult nylon cuff, PVC bladder, PVC inflation bulb with handle, release valve nicked in black.

Supplied in vinyl carrying case with 10 languages user manual and box (GB, FR, IT, ES, PT, DE, GR, RO, HU, Arabic).

Lightweight design, only 135 g