



Product Code	32694
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	Yes
CH-REP	Yes
EAN/UPC	4045396001879
PARAF	951304967

Description

MINIMUS® III sphygmo - twin tube

Minimus® III are high quality, compact and corrosion resistant aneroid sphygmomanometers, featuring a chromium-plated screwed tube connector and ergonomic stainless-steel spoon-shaped grip.

Latex bulb with chromium-plated air release valve.
Precision air release valve with fine adjustment and wearfree.
Linear aluminium graduated-scale Ø 49 mm, up to 300 mmhg.
Maximum error tolerance: ±3 mmhg.

Supplied with adult cuff and vinyl zipper case.

Made in Germany.

Metal connectors for easy and fast changes of cuffs