

# KOBOLD ANEROID SPHYGMOMANOMETER with 3 paediatric cuffs - grey



Product Code	32695
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	178516
CND	C9006
EAN/UPC	4250296155359
UMDNS	16-156
PARAF	987474994

## Description

Elegant and compact, single tube sphygmomanometer for frequent use in clinic and consulting room. Precision manometer, anodised aluminium Ø 48 mm. Sensitively adjustable air reduction valve, high quality chrome finish. Scale: 300 mmHg, 2 mm scale division. New technology with air intake via the manometer capsule in the specially tempered copper-beryllium membrane. Codes 32693, 32695, 32696 are latex free. Made in Germany.

Includes 3 cuffs:  
- Child: 20.5-28 cm  
- Baby: 14-21.5 cm  
- Infant: 10-15 cm