

ERKA SWITCH SIMPLEX ANEROID SPHYGMOMANOMETER - 2 tubes



Product Code	32697
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1147714
CND	C9006
EAN/UPC	4250296160476
UMDNS	16-156
PARAF	987475009

Description

RKA SWITCH SIMPLEX - dial Ø 56 mm - twin tubes grey - 294.408.94

It is manufactured from shock-resistant polycarbonate and has passed the shock test according to AAMI SP9. The sight glass is chemically stable and unbreakable. The brass valve screw works smoothly and is adjustable to exact settings. Measuring range 0-300 mmHg $\pm 1\%$. Black velcro cuff. Supplied with pouch. Made in Germany. 2 year warranty.